

Use and Safety-Instructions

Product: Activator Tubes

These use and safety instructions are valid for Lewa Activator Tubes 315 AR. These tubes are made to work together with a Lewa-Dental Facebow for sagittal movement. Each of our products is functionally tested. However, should a product not function impeccably, despite of our severe controls, we ask you to send it back to us without delay, under mention of the LOT - Number.

1. Insertion

- Keep care that the activator tube is clean, especially free of oil, grease etc.
- Fix the activator tube according the demanded function into the prepared plaster - model by using modelling wax.
- Keep care, that the tube itself is isolated with modelling wax
- Fix the wires and other needed parts at the same manner.
- Spray the acrylic (no Lewa products available) into the model.
Be aware, that all parts, especially the activator tubebody is totally encased of the acrylic.
- Polymerisation of the acrylic according the instructions of the manufacturer.
- Clean the complete plate with hot steam to eliminate the wax
- Check if the arms of the facebow can be put into the tubes.

2. Treatment

Activator Tubes are working together with Facebows. Be aware, that there may be a risk of hurting the patient if the facebows are used without safety device. The treatment should be done under medical supervision. In the case of long treatment, one runs the risk that the parts settle in the mucous membrane. The treating physician must be aware of that and seize preventive measures.

Caution: This product contains Nickel and Chrom! By patients who present allergic reactions (appearance of bubbles, itching), the treatment should be interrupted.

3. Recommendations to the patient

The appliance should be cleaned very carefully. Dirt, lime and tartar can block the function of the appliance. Clean it daily with the toothbrush, one time a week it can be cleaned with hot vinegar water.

In the case of permanently inserted devices, take care of a special mouth hygiene. One should recommend to the patient an intensive clean by means of the toothbrush and dental clean-solution.

