

# NEW! IMPROVED! and IN-STOCK!

## Nitom<sup>2</sup>™ - Laser Welded Locking Facebow

The only locking facebow available has been improved with laser welding, and now comes in it's own zipped headgear pouch. Nitom<sup>2</sup>™ locks into buccal tubes, greatly reducing the risk of injury to patients.

- Easily locked or removed from buccal tubes
- Fits most standard\* tubes .045" (1.13mm)



### Order Form

Please supply me with:

Size	Inner Bow Size	Short Outer Bow	Standard Outer Bow		
		Order No	Qty	Order No	Qty
1	83mm	1200-511	<input type="checkbox"/>	1200-611	<input type="checkbox"/>
2	90mm	1200-512	<input type="checkbox"/>	1200-612	<input type="checkbox"/>
3	97mm	1200-513	<input type="checkbox"/>	1200-613	<input type="checkbox"/>
4	104mm	1200-514	<input type="checkbox"/>	1200-614	<input type="checkbox"/>
5	111mm	1200-515	<input type="checkbox"/>	1200-615	<input type="checkbox"/>

Name: \_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Contact Name: \_\_\_\_\_

Contact Telephone: \_\_\_\_\_

Account Number: \_\_\_\_\_

**Buy 10 at  
only £7.80  
each**

