

# Exhibition Season Offers

**ORTHO**  
**CARE**  
(UK) Limited

## What's New



Ortho-Care (UK) Ltd, 1 Riverside Estate, Saltaire, West Yorkshire, BD17 7DR.  
Tel: +44 (0) 1274 533233 Email: [info@orthocare.co.uk](mailto:info@orthocare.co.uk) . [www.orthocare.co.uk](http://www.orthocare.co.uk)

**Elastic Organiser  
Available Now**





## Aligner & Retainer Cleaning Tablets

| REF                       | List Price |
|---------------------------|------------|
| 9001-239S (On The Go PK4) | 83p        |
| 9001-237 (36 box)         | £4.85      |
| 9001-238 (96 box)         | £9.04      |

**SPECIAL  
OFFER  
55P P/PK**

**BUY 12 BXS  
GET 12 BXS FREE**

Equivalent to  
£2.42

**BUY 6 BXS  
GET 6 BXS FREE**

Equivalent to  
£4.52



Includes a **FREE RETAINER CASE**  
WITH EVERY BOX EVERY TIME! Worth €1.50

## Sonic Cleaner

| REF      | List Price |
|----------|------------|
| 9222-002 | £11.24 ea  |

**BUY 10 at  
£8.37 ea  
OR  
BUY 50 at  
£6.25 ea**



- Universal Cleaning Pod
- Includes 2 F&G Tablets + Batteries

# OrthoCare Relief Wax

| REF       | List Price |
|-----------|------------|
| 860-1110C | £62.90     |



- 5 strips p/pk
- 100 patient pks/box
- Tamper proof seal

**BUY 3 BXS &  
GET 1 BX FREE**

Equivalent to  
£47.18 p/box



Remember, we can personalise your Relief Wax boxes with **your artwork**.

Ask for a quote

# Premium Aligner/Retainer Case with Mirror

| REF       |         | Printed | Non printed |
|-----------|---------|---------|-------------|
| 825-009BK | (black) | £3.98   | £2.88       |
| 825-009WH | (white) | £3.98   | £2.88       |



**BUY 25 GET  
25 FREE**

Printed or  
non printed

1 off set up fee of £26.25 for new artwork

# Orthodontic Mouthwash

| REF                   | List price     |
|-----------------------|----------------|
| 8000-316 (30ml)       | £1.64          |
| 8000-315 (500ml x 12) | £52.42/12 pack |

FREE Orthodontic toothbrush with every 500ml bottle of Fresh & Go Mouthwash! Choose from:

- Travel Brush
- V2 toothbrush
- OrthoSpace toothbrush



**BUY 3 CASES  
GET 1 CASE FREE**  
(Equivalent to £39.32 per case of 12 500ml bottles)  
**+ 48 TOOTHBRUSHES FREE**  
(500ml bottle only)  
**BUY 6 CASES  
GET 3 CASES FREE**  
(Equivalent to £34.95 per case of 12 500ml bottles)  
**+ 108 TOOTHBRUSHES FREE**





# Home-Care Kits



## Fixed Brace

Home-Care Pack

| REF      | List price |
|----------|------------|
| 8000-600 | £8.16      |

**BUY 10 @ £7.70 ea**  
**BUY 25 @ £7.40 ea**



## Deluxe Fixed Brace

Home-Care Pack

| REF      | List price |
|----------|------------|
| 8000-562 | £12.51     |

**BUY 10 @ £11.73 ea**  
**BUY 25 @ £11.03 ea**



## Orthodontic Retainer

Home-Care Pack

| REF          | List price |
|--------------|------------|
| 8000-601-3/4 | £20.07     |

**BUY 10 @ £18.69 ea**  
**BUY 25 @ £17.59 ea**



## Clear Aligner

Home-Care Pack

| REF          | List price |
|--------------|------------|
| 8000-602-3/4 | £22.61     |

**BUY 10 @ £20.90 ea**  
**BUY 25 @ £16.69 ea**



## Mini Aligner

Home-Care Pack

| REF      | List price |
|----------|------------|
| 8000-603 | £4.68      |

**BUY 10 @ £4.36 ea**  
**BUY 25 @ £3.86 ea**

Remember, all our packs can be personalised to suit both **your artwork and your budget.**

*Ask for a quote*





# Bundle Offer



**Aligner & Retainer  
Cleaning Tablets x24**  
Worth £116.40

**Sonic Cleaners x24**  
Worth £269.76

**Printed or Plain  
Retainer Case x24**  
Worth £42.96

**Worth £429.12**

24 | 24 | 24  
**£199.00**

Saving of £230.12

# ToothBrushes

## V2

| REF      | List price |
|----------|------------|
| 8000-542 | £1.66      |



## OrthoSpace

| REF      | List price |
|----------|------------|
| 8000-531 | £1.66      |



## Travel Brush

| REF      | List price |
|----------|------------|
| 8000-521 | £1.66      |



**BUY 36 & GET 12 FREE**

**BUY 72 & GET 36 FREE**

(As little as £1.10ea)

# Aligner Chews

| REF        | List Price                |
|------------|---------------------------|
| 9001-AC-01 | £15.55 (10pks of 2 chews) |



**£8.99 PER PACK  
OF 10 PACKS OF 2**

# Platypus® Orthodontic Toothbrush

| Product                     | REF      | List Price |
|-----------------------------|----------|------------|
| Platypus® Orthodontic Brush | 8000-108 | £3.34      |



**BUY 12 & GET  
12 EXTRA FREE**

# Platypus® Orthodontic Flosser

| Product                       | REF      | List Price |
|-------------------------------|----------|------------|
| Platypus® flossers (pk of 30) | 8000-106 | £5.90      |



**BUY 8 BAGS  
GET 4 EXTRA  
BAGS FREE**

For bulk offers, please ask our Customer Service team.

# BraceSpace Interdental Brushes

| Product  | Size      | REF      | List Price |
|----------|-----------|----------|------------|
| Pink     | 0.4mm     | 8000-260 | £2.77      |
| Orange   | 0.45mm    | 8000-261 | £2.77      |
| Red      | 0.5mm     | 8000-262 | £2.77      |
| Blue     | 0.6mm     | 8000-263 | £2.77      |
| Yellow   | 0.7mm     | 8000-264 | £2.77      |
| Assorted | 1 ea size | 8000-265 | £2.77      |



**BUY 36 & GET 12 FREE**  
Equivalent to £2.08 ea

**BUY 72 & GET  
36 FREE**

Equivalent to £1.85 ea

# Interbrace Brush & Re-fills

Ideal for cleaning around and through braces.

| Product                | REF      | List Price |
|------------------------|----------|------------|
| Interbrace brush       | 8000-255 | £2.92      |
| Refills (12 pks of 10) | 8000-257 | £32.55     |



**BUY 36 BRUSHES & GET 12 FREE**  
Equivalent to £2.19 ea

**BUY 72 BRUSHES & GET 36 FREE**

Equivalent to £1.95 ea

(this offer applies to brushes only not refills)



## Elastic Storage Dispenser

| REF     | List price |
|---------|------------|
| EL-ORG1 | £65.00     |



**SPECIAL  
OFFER  
£45**

For long ligs such  
as Elast-o-loops



Elastics not  
included.

**NEW**

## Spectrum Plus Starter Kit

| REF          | List price |
|--------------|------------|
| 499-400 Kits | £55.07     |



+2 Mosquito Forceps

**SPECIAL  
OFFER  
£45**

| REF       | List price     |
|-----------|----------------|
| 499-Range | £20.64/pk 1000 |



**BUY 3 PKS  
GET 1 PK FREE**  
**BUY 5 PKS  
GET 5 PKS FREE**  
(As little as £10.32 p/pk)

## Elast-o-loops

| REF       | List price     |
|-----------|----------------|
| 466-Range | £39.21/pk 1008 |



**BUY 3 PKS  
GET 1 PK FREE**  
**BUY 5 PKS  
GET 5 PKS FREE**  
(As little as £19.60 p/pk)

## Centric - Elastomeric Chain

| REF        | List Price |
|------------|------------|
| 3077-Range | £16.79/spl |



Closed Link/Continuous



Medium Link/Short Open



Wide Link/Long Open



**BUY 3  
GET 1 FREE**  
**BUY 6  
GET 4 FREE**  
(As little as £10.07 p/pk)

## Dura-Chain Highly Elastic Clear Chain

| REF      | List price |
|----------|------------|
| 3002-849 | £52.49/spl |
| 3004-049 | £52.49/spl |
| 3005-149 | £52.49/spl |



3002-849 - Closed Link/Continuous



3004-049 - Medium Link/Short Open



3005-149 - Wide Link/Long Open



**BUY 3  
GET 1 FREE**  
**BUY 6  
GET 6 FREE**  
(As little as £26.24 p/pk)

# Reliance Products

## Ortho-Flextech™ Dental Retention System

| REF     | Product                      | List Price | OFFER |
|---------|------------------------------|------------|-------|
| R/SSOFT | 76.2cm (30") Stainless steel | £148.83    | £89   |
| R/OFT2  | 76.2cm (30") White Gold      | £320.00    | £299  |



## FlowTain™

| REF   | Product              | List price | OFFER  |
|-------|----------------------|------------|--------|
| R/FT  | Flowable Paste 1.5GM | £42.05     | £26.95 |
| R/FTV | Low Visc Paste 1.5g  | £42.05     | £26.95 |



## GoTo®

Light Cure Bracket Adhesive

| REF   | Product    | List price | OFFER  |
|-------|------------|------------|--------|
| R/GTP | 4g Syringe | £40.79     | £26.95 |



## TKO

| REF    | List price | OFFER  |
|--------|------------|--------|
| R/TKOP | £40.10     | £26.95 |



## Ultra Band-Lok®

Light Cure Band Cement

| REF        | Product                                  | List price | OFFER  |
|------------|------------------------------------------|------------|--------|
| R/NFUBLBPP | 5g Syringe (blue non-fluoride)           | £43.22     | £28.50 |
| R/NFUBLPP  | 5g Syringe (tooth coloured non-fluoride) | £43.22     | £28.50 |



## Assure® PLUS

Bonding Primer

| REF       | Product | List price | OFFER   |
|-----------|---------|------------|---------|
| R/PLUS2ML | 2ml     | £77.14     | £69.42  |
| R/PLUS    | 6ml     | £138.99    | £125.09 |



## Reliance Orthodontic Bonding Kit

| REF       | List price |
|-----------|------------|
| R/BONDKIT | £114.23    |

Kit is worth £143



**SPECIAL OFFER**  
**£99**

## Auxiliaries

Wide range of auxiliaries available

10 p/pk



Scan to see our whole range.

**SPECIAL OFFER**

**£22**

\*Cosmetics not included.



## Gestenco

Freddy Ceramic Hook  
5028-191

Ceramic Bite Bumper  
5028-192

Freddy Ceramic Buttons  
5028-190

**SPECIAL OFFER**

**£45p/pk**

10pk



Scan for more info & higher quality imagery.



| REF      | List price    |
|----------|---------------|
| 5028-190 | £54.50 10p/pk |
| 5028-191 | £54.50 10p/pk |
| 5028-192 | £57.70 10p/pk |



# Molekur™ Molekur™ One S




## NEW PCTG Aligner & Retainer Material

### PERFORMANCES

- Excellent stain-resistant, outstanding properties of anti yellowing and stain resistance.
- It can resist environmental photo degradation and dietary colouring.
- Excellent stress relaxation, great force retention behaviour compared to other similar products on the market.

### PCTG ADVANTAGES

- PCTG has significantly lower breakage rate than PETG.
- PCTG has a lower tensile modulus than PETG, where exhibiting lower/ gentle initial force but better force retention performance than PETG.
- The feature of point 2 also brings better wearing comfort for patients.
- BPA free.



**SPECIAL  
OFFER**  
**96p**  
per blk

SCAN HERE  
for more info



Boxes of 120 or 150 blanks to help reduce packaging waste.

Sold in boxes of 120 blanks  
12 foils of 10

9002-320 - 1mm 125mm Square

9002-321 - 1mm 120mm Round

9002-322 - 1mm 125mm Round

£115.20

Sold in boxes of 150 blanks  
15 foils of 10

9002-317 - 0.76mm 125mm Square

9002-318 - 0.76mm 125mm Round

9002-319 - 0.76mm 120mm Round

£144.00

### REVIEW

"Molekur retainer material is great to work with, forming smoothly without folds or creases and producing a high level of clarity. It is easy to trim using burs or a hot knife, finishes cleanly with Komet brushes, and offers outstanding strength and durability."

Garry Shaw  
Orthodontic Technician

## BUY 2 GET 1 EXTRA FREE



**Distal End Cutter**

| REF    | List price |
|--------|------------|
| 63-101 | £93.45     |



**Distal End Cutter Mini Long Handle**

| REF        | List price |
|------------|------------|
| 63-10105LM | £93.45     |



**Ligature Cutter**

| REF    | List price |
|--------|------------|
| 63-105 | £93.45     |



**Weingart Utility**

| REF    | List price |
|--------|------------|
| 63-201 | £93.45     |



**Angled Weingart Utility**

| REF     | List price |
|---------|------------|
| 63-201A | £93.45     |



**Cinch Back**

| REF    | List price |
|--------|------------|
| 63-100 | £99.75     |



**Bracket Remover**

| REF    | List price |
|--------|------------|
| 63-221 | £93.45     |



**Bracket Remover Angled**

| REF     | List price |
|---------|------------|
| 63-221A | £93.45     |



**Crimpable Hook**

| REF    | List price |
|--------|------------|
| 63-222 | £93.45     |

## SPECIAL OFFERS BELOW ON BONDING INSTRUMENTS



**Bracket Tweezer With Positioner Tip**

| REF    | List price |
|--------|------------|
| 63-928 | £13.60     |

**£8.99**



**Molar Bonding Tweezer Rev Action**

| REF    | List price |
|--------|------------|
| 63-930 | £13.60     |

**£8.99**



**Ortho-Care Mosquito Forcep**

| REF    | List price |
|--------|------------|
| 64-919 | £14.69     |

**£8.99**



**Probe Single Ended No.37 Short**

| REF    | List price |
|--------|------------|
| 60-837 | £6.38      |

**£3.99**



**Mitchells Trimmer**

| REF    | List price |
|--------|------------|
| 60-812 | £16.62     |

**£8.99**



**Task Ligature Director/Tucker**

| REF    | List price |
|--------|------------|
| 60-273 | £32.67     |

**£13.99**

**Mathieu Needle Holder**

| REF      | List price |
|----------|------------|
| 4150-913 | £42.06     |

**£27.50**



## 1/2 Price Archwires

### Euroform®/Centric

| Product (Pks 100)        | Type  | Std price | HALF PRICE | REF   |
|--------------------------|-------|-----------|------------|-------|
| Centric S/E NiTi         | Round | £83.79    | £41.89     | 5823- |
| Centric S/E NiTi         | Rect  | £109.15   | £54.57     | 5823- |
| Centric H/A NiTi         | Round | £87.10    | £43.55     | 5826- |
| Centric H/A NiTi         | Rect  | £131.20   | £65.60     | 5826- |
| Euroform S/S             | Round | £34.18    | £17.09     | 5111- |
| Euroform S/S             | Rect  | £47.41    | £23.70     | 5811- |
| Posted archwires (pk 10) | Rect  | £34.02    | £17.01     | 5322- |

\*Please ask for price on pks 10 Euroform

### Cosmetic Wires

| Cosmetic wires (pk 10) | Type    | Std price | HALF PRICE | REF   |
|------------------------|---------|-----------|------------|-------|
| Euroform NiTi          | Rd/Rect | £48.51    | £24.25     | 5523- |
| Euroform NiTi H/A      | Rd/Rect | £51.93    | £25.96     | 5526- |
| Euroform S/S           | Rd/Rect | £36.38    | £18.19     | 5511- |

### Archwire Holders

| REF       | List price |
|-----------|------------|
| 4811-009  | £64.81     |
| 4811-009S | £64.81     |

**SPECIAL OFFER**  
**£35**



## Titus

Precision Mini Twin Bracket System

| REF      | List price |
|----------|------------|
| 3232-MBT | 72p ea     |



Set =  $\frac{5 - |-5}{5 - |-5}$

| REF        | List price  |
|------------|-------------|
| 3232-Range | 72p/per bkt |



Can be supplied in eco-friendly packaging.

**BUY 50 SETS**  
**GET 25 SETS FREE**  
**U+L 5-5**

**BUY 200 SETS**  
**GET 200 SETS FREE**  
(As little as to 36p per bkt)  
**& receive FREE**  
**Woodpecker**  
(With 200 sets offer)

**Worth**  
**£299.00**



## Signature®

Premium Series Brackets

| REF       | List price |
|-----------|------------|
| 3622-Roth | £2.02 ea   |
| 3722-MBT  | £2.02 ea   |



Set =  $\frac{5 - |-5}{5 - |-5}$

Signature Offer

**BUY 30 SETS**  
**GET 15 SETS FREE**  
**BUY 50 SETS**  
**GET 50 SETS FREE**  
(As little as £1.01)

## Signature® Tubes

1<sup>st</sup>/2<sup>nd</sup> Molar Bondable Tubes

| REF  | 6's | 7's | List price |
|------|-----|-----|------------|
| 351- | 551 | 581 | £2.99 ea   |
| 351- | 552 | 582 | £2.99 ea   |
| 351- | 553 | 583 | £2.99 ea   |
| 351- | 554 | 584 | £2.99 ea   |



Set =  $\frac{6 - |-6}{6 - |-6}$

Signature Tubes Offer  
**BUY 100 @ £1.50 ea**  
**BUY 200 @ £1.29 ea**  
**BUY 400 @ 99p ea**

## Signature® Plus Low Profile

1<sup>st</sup>/2<sup>nd</sup> Molar Bondable Tubes

| REF  | 6's | 7's | List price |
|------|-----|-----|------------|
| 351- | 561 | 591 | £2.99 ea   |
| 351- | 562 | 592 | £2.99 ea   |
| 351- | 563 | 593 | £2.99 ea   |
| 351- | 564 | 594 | £2.99 ea   |



**NEW**



# Glory II Series Self Ligating MBT



| REF            | List Price |
|----------------|------------|
| 3955-000HKS345 | £99.00     |

The Glory II self-ligating metal bracket combines strength, precision, and comfort in a compact design. Made with MIM 17-4PH stainless steel, it offers reliable performance, easy handling, strong bonding, and enhanced patient comfort.



- Smaller overall size
- Larger ligation space
- Tooth number laser marked on the base



## Packaging:

Supplied in a convenient 20 bracket kit

Also available as loose brackets.

**BUY 3 SETS & GET  
1 EXTRA SET FREE**  
As low as £74.25


# Ceramic SL MBT

| REF            | List Price |
|----------------|------------|
| 2022-000HKS345 | £220.00    |

**BUY 3 SETS FOR £555  
& GET 1 FREE SL GATE  
OPENING/CLOSING TOOL**  
Saving of £120

These high-quality self-ligating orthodontic brackets are designed for optimal performance, durability, and aesthetics. Engineered with precision and user convenience in mind, they offer efficient treatment mechanics with reduced friction, while maintaining a smooth and visually appealing finish.

- Monoblock injection ensure smooth surface
- Base design provides strength
- Clear, rounded design



## Packaging:

Supplied in a convenient 20 bracket kit

Also available as loose brackets.

## SL Gate Opening/Closing Tool

| REF        | List price |
|------------|------------|
| 63-SL100-1 | £15.00     |



Open

Close

# Vitrine


| REF        | List Price         |
|------------|--------------------|
| 2055-Range | £119.00 p/full set |

**INTRO OFFER**  
**3 SETS FOR £225**  
 As low as £3.75 P/BKT



Ceramic Brackets made from polycrystalline alumina ensuring clear and transparent. One-body manufacturing injection process solves the problem of brittleness and breakage and guarantees the precise slot.

- High breakage resistance and torque strength
- Good biocompatibility
- Colour stability and no discolouration



**Packaging:**  
 Supplied in a convenient 20 bracket kit

*Also available as loose brackets.*

## Purity® Mini

MBT only  
 Sapphire Crystal

| REF         | List price |
|-------------|------------|
| 2068-000U/L | £286.40    |

## Purity® Plus

MBT only  
 Sapphire Crystal

| REF         | List price |
|-------------|------------|
| 2067-000U/L | £286.40    |

## Purity®


Roth/MBT  
 Sapphire Crystal

| REF         | List price |
|-------------|------------|
| 2077-000U/L | £286.40    |

## Prism™

Roth/MBT  
 Co-Polymer Resin

| REF         | List price |
|-------------|------------|
| 2088-000HKS | £286.40    |



**BUY 5 SETS**  
 @ £7.96 per bkt

**BUY 10 SETS**  
 @ £4.95 per bkt  
 £99 per set

**BUY 5 SETS**  
 @ £7.96 per bkt

**BUY 10 SETS**  
 @ £4.95 per bkt  
 £99 per set

**BUY 5 SETS**  
 @ £7.96 per bkt

**BUY 10 SETS**  
 @ £4.95 per bkt  
 £99 per set

**BUY 5 SETS**  
 @ £7.96 per bkt

**BUY 10 SETS**  
 @ £4.95 per bkt  
 £99 per set

All Bracket  
 Sets = 5-|-5  
 5-|-5



# Dentatus IPR Kit

## 5 Piece File Kit + Hand Piece

| REF       | List price |
|-----------|------------|
| PDX-P/SK2 | £733.95    |

**SPECIAL OFFER**  
**£599**

## 5 Piece File Kit

| REF       | List price |
|-----------|------------|
| LTA-P/SK2 | £261.45    |

**SPECIAL OFFER**  
**£220**

## IPR Tips

| REF      | List price |
|----------|------------|
| LTA-P1-5 | £71.40     |

**SPECIAL OFFER**  
**£65**  
per pk



0.1mm  
2 Pack  
LTA-P1/2



0.2mm  
2 Pack  
LTA-P2/2



0.3mm  
2 Pack  
LTA-P3/2



0.4mm  
2 Pack  
LTA-P4/2



0.5mm  
2 Pack  
LTA-P5/2



## Review

"The Dentatus IPR handpiece is a great bit of kit that makes IPR simple and stress free. The dedicated handpiece provides good control whilst removing tooth tissue and the graded colour coded tips are easy to work through. I find that the calibrated tips allow accurate IPR to be undertaken to the desired prescription without needing to check repeatedly with a gauge."

Liam Monaghan | North West Orthodontics

## IPR Kit

| REF       | List price |
|-----------|------------|
| 9001-5KIT | £353.78    |



**SPECIAL OFFER**  
**£260**  
per kit

## Incremental Thickness Gauge

| REF      | List price |
|----------|------------|
| 9001-240 | £55.83     |



Set of 6 keys

**SPECIAL OFFER**  
**£25**

## Horico Discs

| Product                    | REF            | List price |
|----------------------------|----------------|------------|
| D/S 0.15 DIAM DISC RA      | 9001-W355F-190 | £24.26 ea  |
| S/S 0.10 DIAM DISC RA      | 9001-W356F-190 | £24.26 ea  |
| D/S 0.15 DIAMOND DISC RA   | 9001-W401F-190 | £23.53 ea  |
| D/S 0.17 DIAM DISC 14mm RA | 9001-W934-140  | £28.89 ea  |
| D/S 0.17 DIAM DISC 19MM RA | 9001-W934-190  | £28.89 ea  |
| D/S 0.17 DIAM DISC 22MM RA | 9001-W934-220  | £28.89 ea  |

**£16.99** Applies to all  
Horico discs  
on this page

Double sided



9001-W355F-190

Single sided



9001-W356F-190

Double sided



9001-W934-140

Double sided



9001-W934-190

Double sided



9001-W934-220

Double sided



9001-W401F-190

| Product                                  | REF               | List price        | OFFER         |
|------------------------------------------|-------------------|-------------------|---------------|
| POST ARS X CUT BUR 1.0MM HEAD 5.9 LENGTH | 9001-699C         | £71.61 5pk        | <b>£54.56</b> |
| <b>DIAMOND BUR FOR ARS ANTERIOR</b>      | <b>9001-55000</b> | <b>£55.07 5pk</b> | <b>£36.72</b> |
| SAFE TIPPED DIAMOND BURS FOR STAR        | 9001-5MDST        | £78.83 5pk        | <b>£60.06</b> |
| SAFE TIPPED DIAMOND BURS FOR STARS       | 9001-5FDST        | £78.83 5pk        | <b>£60.06</b> |
| X CUT TAPERED FISSURE ARS CARBIDE BUR    | 9001-699LC        | £55.07 5pk        | <b>£41.96</b> |
| FINE GRIT DIAMOND ARS BUR                | 9001-514FD        | £49.17 5pk        | <b>£37.46</b> |



= Best Seller

## IPR Bur Kit

| REF      | List price |
|----------|------------|
| 9001-510 | £146.95    |



## IPR Disc Kit

| REF      | List price |
|----------|------------|
| 9001-511 | £146.95    |



**£139.95**

Applies to both kits on this page

For more info on the benefits and kit contents **SCAN HERE**



## Komet Diamond Strips

| REF              | Product                                                | Qty  | List price |
|------------------|--------------------------------------------------------|------|------------|
| DS25             | Komet Vision-flex 2.5mm Medium Single Sided 45m Blue   | PK10 | £93.71     |
| DS25F            | Komet Vision-flex 2.5mm Fine Single Sided 30m Red      | PK10 | £93.71     |
| <b>NEW</b> DS25A | Komet Vision-flex 2.5mm Assorted Single Sided          | PK15 | £132.30    |
| DS37D            | Komet Vision-flex S3.75mm Medium Double Sided 45m      | PK10 | £93.71     |
| <b>★</b> DS37    | Komet Vision-flex S3.75mm Medium Single Sided 45m Blue | PK10 | £93.71     |
| DS37F            | Komet Vision-flex S3.75mm Fine Single Sided 30m Red    | PK10 | £93.71     |
| <b>NEW</b> DS37A | Komet Vision-flex S3.75mm Assorted Single Sided        | PK15 | £132.30    |

### OFFER

**£65**

**£65**

**£99**

**£65**

**£65**

**£65**

**£99**

## Komet Diamond Strips - Perforated/Honeycomb

| REF           | Product                                                        | Qty  | List price |
|---------------|----------------------------------------------------------------|------|------------|
| WS25          | Komet Vision-flex Narrow 2.5mm Medium Grit Blue Single Sided   | PK10 | £104.74    |
| WS25A         | Komet Vision-flex Narrow 2.5mm Assorted Single Sided           | PK15 | £132.30    |
| WS25EF        | Komet Vision-flex Narrow 2.5mm X-Fine Grit Yellow              | PK10 | £104.74    |
| WS25F         | Komet Vision-flex Narrow 2.5mm Fine Grit Red Single Sided      | PK10 | £104.74    |
| <b>★</b> WS37 | Komet Vision-flex Single Sided Strips Wide 3.7mm 0.13mm Blue   | PK10 | £104.74    |
| WS37A         | Komet Vision-flex Single Sided Strips Wide 3.7mm Assorted      | PK15 | £132.30    |
| WS37EF        | Komet Vision-flex Single Sided Strips Wide 3.7mm 0.08mm Yellow | PK10 | £104.74    |
| WS37F         | Komet Vision-flex Single Sided Strips Wide 3.7mm 0.10mm Red    | PK10 | £104.74    |

**£79**

**£99**

**£79**

**£79**

**£79**

**£99**

**£79**

**£79**

## Edenta StripFIX Kit

| REF          | List price |
|--------------|------------|
| 9001-SF-SET1 | £75.17     |



**SPECIAL OFFER**  
**£48**  
per kit

## Edenta Diamond Perforated Strips

| Product (pack of 10)            | REF    | List price |
|---------------------------------|--------|------------|
| <b>★</b> Blue wide s/s med grit | PXDS3  | £57.77     |
| Yellow wide s/s x-fine grit     | PCXDS3 | £57.77     |
| Red wide s/s fine grit          | PFXDS3 | £57.77     |
| Red narrow s/s fine grit        | PFXDS2 | £57.77     |
| Yellow narrow s/s x-fine        | PCXDS2 | £57.77     |
| Blue narrow s/s med grit        | PXDS2  | £57.77     |

## Edenta Diamond Finishing Strips

| Product (pack of 5)                           | REF   | List price |
|-----------------------------------------------|-------|------------|
| Yellow double sided X-fine grit 0.10mm        | KDSC1 | £50.35     |
| Red double sided fine grit 0.20mm             | KDSF2 | £50.35     |
| Red double sided fine grit 0.30mm             | KDSF3 | £50.35     |
| <b>★</b> Blue double sided medium grit 0.30mm | KDSM3 | £50.35     |
| Blue double sided medium grit 0.40mm          | KDSM4 | £50.35     |



**£38**  
perforated strips

**£35**  
diamond strips

**★** = Best Seller



# SAFEJAWZ®

**ORTHO SERIES**  
BRACES MOUTHGUARD

| Product       | REF      | List price |
|---------------|----------|------------|
| Ice blue      | 5100-021 | £18.22     |
| Ice pink      | 5100-022 | £18.22     |
| Crystal clear | 5100-024 | £18.22     |
| Red           | 5100-025 | £18.22     |
| Black         | 5100-026 | £18.22     |

**BUY 10 & GET 4 FREE**  
**BUY 20 & GET 10 FREE**  
Offer applies to all SAFEJAWZ products  
(As little as £12.15)



**BRACEHARBOUR™**



**FLUIDFIT™**

**Fluid Fit™** technology means that SAFEJAWZ® mouthguards feature a super slim profile and will clip on to the teeth for a comfortable fit so that breathing and speech will not be a problem.



**JAWSECURE™**

The Ortho Series Bracket Mouthguard has all the impressive features of the standard SAFEJAWZ®, but with special design features that enable bracket wearers to have a great fitting mouthguard, PROVIDING UPPER AND LOWER PROTECTION FOR TEETH, GUMS AND LIPS ALL THROUGH TREATMENT.

**Jaw Secure™** offers a pre-contoured base which can also assist with stabilising the jaw and protecting the bottom set of teeth.



**LOWERBLOK™**



**REMODEL TECH™**

**Remodel Tech™** allows the wearer to remodel multiple times as and when teeth move with orthodontic treatment ensuring the best fit at all times. Remodel Tech™ is a trademark exclusive to SAFEJAWZ®.

## INTRODUCING THE NEW MARVEL RANGE



| Product         | REF      | List price |
|-----------------|----------|------------|
| Spider-Man      | 5100-027 | £18.22     |
| Captain America | 5100-028 | £18.22     |
| Venom           | 5100-029 | £18.22     |

**BUY 10 & GET 4 FREE**  
**BUY 20 & GET 10 FREE**  
Offer applies to all SAFEJAWZ products  
(As little as £12.15)